

SMART SENSE

BY DIGI[®]

Retail Pharmacies

Critical Asset Protection for Automated Compliance and Patient Safety

SmartSense is your trusted asset protection platform that continuously monitors temperature and humidity levels to ensure customer safety, regulatory compliance, and operational efficiency. More than simply recording data and producing reports, SmartSense equips frontline workers with critical alerts, guidance for corrective actions, and confirmation of execution. The centralized platform helps you maintain medication efficacy while reducing loss with cellular-based temperature monitoring that surpasses SBOP, CDC, and VFC compliance standards.

A Complete Retail Pharmacy Sensing-as-a-Service Platform



Document and report on asset performance using real-time enterprise-wide monitoring and corrective action management.



Avoid loss of valuable medication with real-time excursion alerts, location services, and notification hierarchies for key stakeholders.



Exceed regulatory requirements for VFC storage and handling and twice daily logging requirements.



Quickly deploy in retail environments using wire-free and cellular-based sensors and gateways that don't rely on store networks.



Ensure compliance with NIST certified sensors, documented proof of temperature and humidity levels, and simple NIST recalibration.



Protect products while in transit with track and trace plus condition monitoring that detects unsafe conditions and tampering.

Protecting Patients and Billions in Inventory

"We have several billion dollars' worth of cold chain product that we dispense in-house in our pharmacies annually. So we needed a top of the line solution to make sure that we had really tight monitoring of that inventory to promote the highest levels of patient safety. We use SmartSense in every single refrigerator or freezer that we have in our pharmacies."



Molly Gombos
Lead Director,
Retail Quality Standards



Simple NIST Recalibration Process

With SmartSense, our 3-step NIST recalibration process is quick and easy:

1. A new calibrated NIST probe is sent to you.
2. Simply replace the old NIST probe. It's a painless twist on recalibration.
3. The new NIST probe automatically uploads the new calibration certificates to SmartSense.

The process takes seconds, resulting in zero interruption to temperature readings. Never again return a sensor to the manufacturer for recalibration.



SmartSense NIST Probes

Vaccine Refrigerator



For more information about SmartSense visit www.smartsense.co or call +1 (866) 806-2653

SMART SENSE
BY DIGI